



Ref.: 94230
SHORT 51cm

Ref.: 94240
LONG 69cm

F ... C .. S I

KNEE ORTHESIS WITH FLEXION-EXTENSION JOINT

SIZES ONE SIZE



Polycentric articulations available in 2 lengths

Characteristics:

Knee orthosis that comprises two foam corsets, two Velcro straps and polycentric articulations adjustable at 0-15-30-60 and 90° to limit the flexion/extension.

Indications:

- For unstable or injured knees.
- Post-operative and post-traumatic rehabilitation.
- For relative immobilisation of the articulation in situations that require perfect control.
- To obtain a range of control of the mobility of the operated knee (crossed ligament surgery).

Colour: ■ Black.

Models:

- 94230: Short.
- 94240: Long.

Kneortose post- op Kort utgave har vi på lager

Art : 94230

Størrelse: one size

Fra produsenten: Orliman



blatchford
ORTOPEDI AS

Tlf: 75 50 74 20

Fax: 75 50 74 21

Org.nr.: 914 120 349 MVA

Besøksadresse:

Gidsken Jakobsens vei 32
8008 Bodø

Postadresse:

Postboks 1114
8001 Bodø

Epost: post.bodo@blatchford.no

Web: www.blatchford.no